


Originlab Installation Guide – Windows

Download OriginLab_v9.4.iso from ISTF website under Software category.


Extract the downloaded iso file on your local system.

Navigate to the location where iso has been extracted.

C:\Users\SP\Desktop\OriginLab_v9.4



Just right click on Setup file and Select Run as administrator.




Installation is in Process.




Click on Next.



Click on Yes.



Select Install Product (requires serial Number) Click on Next.



Enter user name: **IIT**

Enter Company name: **IITGN**

Enter user name: **GF3S5-3089-7907446**

click on Next.

Origin Setup

Customer Information
Please enter your information.

Please enter your name, the name of the company for which you work and the product serial number.

User Name:
IIT

Company Name:
IITGN

Serial Number (enter full Serial Number using format xxxxx-xxxx-xxxxxxx):
GF3S5-3089-7907446

InstallShield

< Back Next > Cancel

Click on Yes.

Registration Confirmation

You have provided the following registration information:

Name:	IIT
Company:	IITGN
Serial Number:	GF3S5-3089-7907446


Is this registration information correct?


Yes No

Select you OS Bit and click on Next.




Click on Next.






Click on Next






Click on Next.




Click on Next.




Installation in process it will take 15 to 20 minutes(Approx.)




Click on Finish.



When you start first time you must run as administrator.



Click on Yes.



Enter Server IP: 10.0.137.114

Enter TCP/IP Port: 28000

Click on Update.

